FILTER FUNNEL STANDS 1300 FILTER FUNNEL STAND

The Model 1300 Filter Funnel Stand, the Model 1390 Funnels, and Model 1397 Sample Bottles provide a simple, efficient system for obtaining soil extracts. By connecting a vacuum source to the Funnel Stand, the sample bottles are evacuated, which draws the extract from the soil and into the bottle.

1300, Filter Funnel Stand (shown with 1390V1 funnels and 1395 1 oz. bottles)

Product No.	Description	Weight
1300	FILTER FUNNEL STAND	0.60 kgs
	Holds up to five 80 ml or three 375 ml Buchner funnels	
Accessories		
Product No.	Description	Weight
1390V1	80 ML BUCHNER FUNNEL, Polypropylene	$0.02 \ \mathrm{kgs}$
	Accepts 5.5 cm diameter filter paper	
1390V5	375 ML BUCHNER FUNNEL, Polypropylene	0.08 kgs
	Accepts 9 cm diameter filter paper	
1395	1 OZ. BOTTLES & CAPS, one dozen	$0.59 \ \mathrm{kgs}$
1396	2 OZ. BOTTLES & CAPS, one dozen	$0.86 \ \mathrm{kgs}$
1397	4 OZ. BOTTLES & CAPS, one dozen	$1.36 \ \mathrm{kgs}$

REPLACEMENT PARTS, FILTER FUNNEL STANDS on page 45